

GENESIS

THE WORLD'S MOST ADVANCED COUNTERTOP WATER DISTILLER



**A TRUE BREAKTHROUGH IN WATER PURIFICATION TECHNOLOGY:
ELIMINATES WATER CONTAMINANTS**

Flash vaporizes tap water, killing and removing over 99% of most contaminants.


GREAT TASTING WATER

Patented distillation technology combined with advanced carbon post-filtration ensures consistently high-quality, fresh tasting, purified water virtually free of all contaminants.

PORTABLE AND EASY TO USE

Just plug it in. No plumbing, wiring or installation required.

*Product open
with boiling tray visible*



GENESISTM

THE WORLD'S MOST ADVANCED COUNTERTOP WATER DISTILLER

PERFORMANCE DATA SHEET- CONTAMINANT REMOVAL DATA

IMPORTANT NOTICE: Read this Performance Data Sheet and compare the capabilities of this unit with your actual water treatment needs. It is recommended that, before purchasing a water treatment unit, you have your water supply tested to determine your actual water treatment need.

GENESIS WATER DISTILLER AND PURIFIER-MODEL 6000

REMOVAL				ACTION LEVEL/ MAXIMUM CONTAMINANT LEVEL
Substance	Influent Water mg/L	Distiller Effluent Avg./Max.	Efficiency Avg./Min.	
Total Dissolved Solids	1050	4.0/5.4 mg/L	99.6%/99.5%	500 mg/L
Arsenic ²	00.30	0.003 mg/L	>99%	00.05 mg/L
Barium ²	10.0	0.1 mg/L	>99%	02.0 mg/L
Cadmium ²	00.03	0.0003mg/L	>99%	00.005 mg/L
Chromium ²	00.3	0.002mg/L	>99%	00.1 mg/L
Copper ²	04.0	0.04mg/L	>99%	01.3 mg/L
Selenium ²	00.10	0.00009mg/L	>99%	00.05 mg/L
Lead ¹	00.17	0.001/0.005	98.8%/97%	00.015 mg/L
Nitrate ¹	33	0.12/0.03	>99%/99%	10 mg/L

1. Performance testing conducted and certified by Spectrum Labs, Inc., New Brighton, MN based on test protocol published by the National Sanitation Foundation (NSF) with specific reference to NSF Standard Number 62, "Drinking Water Distillation Systems", Appendix B, "Contaminant Reduction Test Methods".

2. Based on the study "Evaluation of Total Dissolved Solids as a Surrogate Parameter for the Reduction of Inorganic Contaminants by Distillation Systems," conducted for the Water Quality Association by NSF International. TDS may be used as a surrogate for verifying the reduction of arsenic, barium, cadmium, chromium, copper, lead nitrate, and selenium to equal or below MCL when tested in accordance with Annex B. Refer to Annex F for rationale and supporting data.

US Patent No. 4,818,344 European Patent No. EP 0320804B1 E.P.A. Reg. No. 67629-11K-001
Other Worldwide Patents Issued and Pending
©1994 Genesis. An Environmental Corp. All rights reserved.

FOR SALES INFORMATION CONTACT:

GENESIS 6000 ADVANCED DISTILLER AND PURIFIER- INFORMATION SPECIFICATIONS:

- Requires no assembly.
- Includes two 3 liter/ quart glass containers.
- Stainless steel heating element and condenser coil.
- Fully automatic shut-off.
- 25cm/10" high-rise steam lift.
- Cool distillate - ambient temperature.
- Water pre-heated and vented to remove volatile organic chemicals.
- Optional post-carbon block filtration.
- Removable (slide-out) dishwasher-safe stainless steel lined boiling tray.
- Limited one year warranty.
- Lightweight: 7 kg/ 15 lb.
- Compact Size:
51.5cm W x 33.5cm D x 39cm H
(20" W x 13" D x 15" H).
- Fast Distilled Water
Production: Approximately
5.5 ltr./ 1.5 gallons in 5 hours.
27 ltr./ 7 gallons in 24 hours.
- Certified to meet CSA, UL and CE standards.
120V 60Hz./950 watts; 220V-
240V 50Hz./840-1000 watts.
- Starts producing distilled water in 5 minutes.

